

# Integra™ Jarit®

Slide Lock Laparoscopic Instruments

Limits uncertainty with  
smooth action, precise grip  
and ultimate control



INTEGRA™  
LIMIT UNCERTAINTY

### Grasping Forceps

Straight, serrated jaws with proximal recess, double action

#### 610-100

ROTO-LOK™, slide lock, insulated shaft, 5mm, 32cm

#### 605-100

Slide lock, insulated shaft, 5mm, 32cm



### Grabber Forceps

Long, fenestrated jaws with interlocking wave-like serrations for strong hold, double action

#### 610-104

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm



### Tenaculum Forceps

Single-toothed long jaws designed for grasp of thick tissue, double action

#### 605-152

Slide lock, insulated shaft, 5mm, 32cm



### Fundus Grasping Forceps

Long, beveled-edge jaws permit a large purchase on tissue and maintain a strong hold

#### 610-107 - Single action

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm



#### 610-107DA - Double action

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm



### Dual-Cup Grasping Forceps

Double cup design increases surface area contact which improves hold on tissue, double action

#### 610-108

ROTO-LOK, slide lock, insulated shaft, 5mm, 32cm



smooth action  
precise grip  
ultimate control



## DeBakey Forceps

Fenestrated 40mm long jaws with DeBakey-type teeth, double action

**610-140**

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm

**610-140L**

ROTO-LOK, slide lock, insulated shaft, 5mm, 45cm



## Babcock Forceps

Traditional design for grasping delicate tissue, double action

**610-147**

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm



## Swanström¹ Grasping Forceps

Unique design of 1x2 DeBakey-type teeth with triangular fenestration provides an effective hold with atraumatic grasp, double action

**610-148**

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm



## Strong Hold Babcock Forceps

Traditional design for atraumatic grasping, double action

**610-149**

ROTO-LOCK, slide lock, insulated shaft, 5mm, 37cm



## Hunter² Bowel Grasper

Fenestrated, 20mm jaws with DeBakey-type teeth for atraumatic manipulation of delicate tissue, e.g., bowel, stomach, double action

**610-163**

ROTO-LOK, slide lock, insulated shaft, 5mm, 32cm

**610-164**

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm

**610-165**

ROTO-LOK, slide lock, insulated shaft, 5mm, 45cm



**610-166**

ROTO-LOK, slide lock, insulated shaft, 5mm, 37cm, 45mm jaw

1. Developed in cooperation with Lee L. Swanström, M.D., F.A.C.S., Legacy Portland Hospital, Portland, Oregon.

2. Developed in cooperation with John Hunter, M.D., FACS, Professor and Chairman, Department of Surgery, Oregon Health & Science University, Portland, Oregon.

### Full Catalog Listing

| Catalog Number | Description                                                                               |
|----------------|-------------------------------------------------------------------------------------------|
| 605-100        | Grasping Forceps Double action, insulated shaft, 5MM, 32CM                                |
| 610-100        | Grasping Forceps ROTO-LOK, double action, insulated shaft, 5MM, 32CM                      |
| 610-107        | Fundus Grasping Forceps ROTO-LOK, single action, insulated shaft, 5MM, 37CM               |
| 610-107DA      | Fundus Grasping Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM               |
| 610-163        | Hunter Bowel Grasper ROTO-LOK, double action, insulated shaft, 5MM, 32CM                  |
| 610-164        | Hunter Bowel Grasper ROTO-LOK, double action, insulated shaft, 5MM, 37CM                  |
| 610-165        | Hunter Bowel Grasper ROTO-LOK, double action, insulated shaft, 5MM, 45CM                  |
| 610-166        | Hunter Bowel Grasper ROTO-LOK, double action, insulated shaft, 5MM, 37CM, 45mm jaw length |
| 610-140        | DeBakey Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM                       |
| 610-140L       | DeBakey Forceps ROTO-LOK, double action, insulated shaft, 5MM, 45CM                       |
| 605-152        | Tenaculum Forceps Double action, insulated shaft, 5MM, 32CM                               |
| 610-104        | Grabber Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM                       |
| 610-108        | Dual-Cup Grasping Forceps ROTO-LOK, double action, insulated shaft, 5MM, 32CM             |
| 610-147        | Babcock Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM                       |
| 610-148        | Swamström Grasping Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM            |
| 610-149        | Strong Hold Babcock Forceps ROTO-LOK, double action, insulated shaft, 5MM, 37CM           |

For more information or to place an order, please contact:

Integra ■ 589 Davies Drive, York, PA 17402

877-468-5572 USA ■ 914-789-7094 outside USA

609-750-4259 fax ■ 914-592-8056 outside USA fax

[integralife.com](http://integralife.com)