




ACF/MLD RETRACTOR SYSTEM



Multiple Options For Multiple Needs



ACF/MLD RETRACTOR SYSTEM

VERSATILITY...

STRENGTH...

EASE OF USE.

With multiple blade and retractor options, the Redmond R&B™ ACF/MLD Retractor System is used for single-level and multi-level procedures in cervical fusion and micro lumbar discectomy cases. The system remains secure during the entire procedure and is designed for easy identification, insertion, and retraction. The ACF/MLD Retractor System meets the numerous needs of today's spine surgeon.

Single-Level Configuration



Multi-Level Configuration

ACF/MLD RETRACTOR SYSTEM

- Dual blade attachment points allow for expanded retraction to perform single-level and multi-level surgeries.
- Transverse and Longitudinal Retractors for 4-sided retraction.
- Multiple hinge points to match anatomical variations.
- Specially-designed curve in retractor arms for easier blade attachment.
- Blades attach to retractor arm along entire width for a more secure fit.



Longitudinal Retractor



Transverse Retractor

- Spring-loaded handle for easy insertion and release.



Retractor Handle

- Color-coded depth gauge to eliminate guesswork.



Depth Gauge

- Color-coded blades and hooks for easy identification.
- 25mm to 75mm lengths for multiple anatomies and procedures.
- Aluminum construction for high strength with radiolucency.




Short-teeth Blades

- Multiple blade styles:
 - 15mm and 20mm widths
 - Hooks, blunt blades, short teeth blades and long teeth blades



Hook and Blade Assortment

ACF/MLD RETRACTOR SYSTEM



Catalog #	Description	Catalog #	Description
RB2400	Complete set including 2 retractor handles and sterilization case	RB2405	Blunt blade set, 12 pair
RB2408	Transverse retractor body, double level	RB2406	Long teeth blade set, 12 pair
RB2409	Longitudinal retractor body, double level	RB2407	Hook set, 6 hooks
RB2403A	Retractor handle (2 per set)	RB2422	ACF/MLD Depth Gauge (color-coded)
RB2404	Short teeth blade set, 12 pair	RBT2400	ACF/MLD Sterilization case



J. Jamner Surgical Instruments
 9 Skyline Drive • Hawthorne, NY 10532
 (914) 592-9050 • (800) 431-1123 • Fax: (914) 592-8056

R&B Redmond is a trademark of Integra LifeSciences Corporation.
 J. Jamner Surgical Instruments Inc., is a subsidiary of Integra LifeSciences Corporation.
 ©2004 J. Jamner Surgical Instruments, Inc. All Rights Reserved.

NS642 (9/04)
 Printed in USA
 2k