

Miltex® Surgical Instrument Stain Remover

Powder formula removes stains, tarnish and restores luster to surgical instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, etc. Ideal for use on stainless steel, chrome, nickel, brass, copper and silver.

- Effective on instruments, autoclaves, ultrasonic cleaners or any stainless steel, tile or plated surface
- Will not harm instrument surfaces
- Easy, fast, convenient to use
- Economical - a little goes a long way




REF	Description
3-740	Surgical Instrument Stain Remover, 3 oz. (85 grams), (case of 12)*

*(Contact Distributor for availability of individual units)

Miltex® Nylon & Stainless Steel Brushes

- Nylon bristles for effective cleaning of instrument box locks
- Stainless steel bristles for cleaning serrations
- Autoclavable



REF	Description
3-1000	Instrument Cleaning Brushes 7-1/4" (18.4cm), Nylon Bristles (Set of 3)
3-1001	Instrument Cleaning Brushes 7-1/4" (18.4cm), Stainless Steel Bristles (Set of 3)
21-379	Stainless Wire Brush, 6-1/4" (15.9 cm)

Miltex® Instrument Care Kit System

A great way to try various Miltex® products!

REF	Description
3-800	Spray Lube (3-700), Surgical Instrument Cleaner (3-720), Surgical Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon Brush and Stainless Steel Brush



Miltex® Instrument Care Training

25-minute Instructional DVD assists practitioners and staff to properly care for instruments and to keep them in top working condition. DVD offers step-by-step guidance on handling, cleaning, rinsing, inspecting, lubricating and sterilizing surgical instruments.

Step-by-step poster provides easy instructions to clean and prepare instruments for sterilization.

Try a variety of our popular cleaning products with a H.I.P. sample box.

REF	Description
3-810	Miltex® Instrument Care DVD
SURGHIPPOST	Miltex® H.I.P. Poster
SURGHIPBOX	Miltex® H.I.P. Sample Box



SURGHIPPOST



3-810



SURGHIPBOX

www.miltex.com
toll-free phone 866-854-8300
toll-free fax 866-854-8400
phone 717-840-9335
fax 717-840-9347

Please call Customer Service today for additional information, an in-office demonstration, to place an order, or to receive other Miltex catalogs.

Copyright © 2010 Integra LifeSciences Corporation All Rights Reserved. Miltex and EZ-Zyme are registered trademarks and Miltex Redefining Excellence and Tip-It are trademarks of Miltex, Inc. a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra Life Sciences Corporation.

0410/10/H/9

MPM0506

INST-CARE-BROC



+H834INSTCAREBROC07



INSTRUMENT CARE & CLEANING PRODUCTS

Get H.I.P. with the MilteX® Healthy Instrument Program!

The MilteX® care and cleaning product portfolio is designed to provide optimum cleaning performance and will help to prolong the service life of valuable surgical and dental instrumentation.

MilteX® Instrument Prep Enzyme Foam

Ready-to-use foaming spray for pre-cleaning of soiled instruments and scopes (rigid or flexible). Ideal for applications where instrument cleaning is not immediately available. Simply spray on instruments until ready for decontamination and cleaning.



REF	Description
3-760	Instrument Prep Enzyme Foam, 24 oz. (0.71 liter) spray bottles (case of 12)*

*(Contact Distributor for availability of individual units)

MilteX® EZ-Zyme® All Purpose Enzyme Cleaner

Ideal for cleaning hard to remove proteinaceous (mucus & tissue) and lipid (fatty) materials from the surfaces of surgical and dental instruments, linen, water submersible equipment, or where mechanical scrubbing is impractical.

- Easy to Use
- Neutral pH
- Multi-use as an instrument presoak, ultrasonic cleaner, laundry presoak and evacuation system cleaner
- Space-saving or bulk packaging options



REF	Description
3-750	EZ-Zyme® All Purpose Enzyme Cleaner, 3/4 oz. (22 ml) packets (dispenser box of 32). Contents of 1-box makes 32 gallons of solution
3-755	EZ-Zyme® All Purpose Enzyme Cleaner, 1-gallon bottle (3.79 liter). Contents of 1-gallon container makes 170 gallons of solution

MilteX® Surgical Instrument Cleaner & Lubricant

A dual-purpose instrument cleaner and lubricant. Chemically engineered for use as a presoak or for manual cleaning, washer/decontaminators, ultrasonic cleaners and cart washers.



- Cleaner & Brighter Instruments
- Removes Stains, Blood and Mineral Deposits
- Biodegradable
- Concentrate, Lower Cost Per Use

REF	Description
3-710	Surgical Instrument Cleaner and Lubricant 1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

MilteX® Surgical Instrument Cleaner

Super concentrate formula for manual and ultrasonic cleaning or soaking of surgical instruments, accessories and glassware.

- Neutral pH
- Biodegradable
- Rinses free, leaves no residue
- Phosphate-free
- Concentrated



REF	Description
3-720	Surgical Instrument Cleaner, 8 fl. oz. (0.24 liter) bottles (case of 12). Contents of one 8 oz. bottle makes 32 gallons of solution*
3-725	Surgical Instrument Cleaner, 1-gallon bottle (3.79 liter). Contents of 1-gallon container makes 512 gallons of solution
3-726	Surgical Instrument Cleaner, 1/4 fl. oz. (7.4 ml) packets (dispenser box of 80). Contents of 1-box makes 80 gallons of solution

*(Contact Distributor for availability of individual units)

MilteX® Spray Lube (Surgical Instrument Lubricant)

Pre-mixed, non-silicone and anti-corrosive formula in a ready to use 8 oz. spray pump bottle. Steam permeable and water soluble. Helps prevent spotting, staining and rusting. Lubricates instrument box locks, screw locks, scissors blades and any other metal on metal moving part.



- Ideal for spot treatment

NOTE: This product differs from the MilteX® Instrument Lubricant (REF 3-705) and is not intended for instrument immersion.

REF	Description
3-700	Spray Lube, 8 fl. oz. (0.24 liter) bottles (case of 12)*

*(Contact Distributor for availability of individual units)

MilteX® Instrument Lubricant

A water emulsible soaking agent works to extend the life of surgical instruments. Concentrated non-silicone lubricant and rust inhibitor. Use daily to keep moving parts operating smoothly.



- Formulated to extend the working life of surgical instruments
- Helps Prevent Corrosion
- Formulated with gentle, non-aggressive ingredients
- Frees Up Box Locks - Lubricates moving parts of surgical instruments

REF	Description
3-705	Instrument Lubricant, 1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

MilteX® Tip-It™ Instrument Protectors

MilteX®Tip-It™ Instrument Protectors guard instrument tips and edges from damage during handling and sterilization.

- PVC medical grade plastic
- Latex-free material
- Large holes (vented styles) to reduce vacuum and moisture build up for easy removal from instrument
- Steam and gas sterilizable
- 50 per pack



REF	REF	REF	REF	Dimension	Size Code	Color
Clear Vented	Tinted Vented	Solid Vented	Solid Non-Vented			
3-2501C	3-2501T	3-2501V	3-2501	1/16" X 3/4" (1.6mm X 19mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" X 3/4" (2mm X 19mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" X 3/4" (2.8mm X 19mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" X 1" (3.2mm X 25.4mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" X 1" (4.8mm X 25.4mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" X 3/4" (10.3mm X 19mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" X 3/8" X 1" (1.6mm X 9.5mm X 25.4mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" X 5/8" X 1" (2mm X 16mm X 25.4mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" X 1" X 1" (3.2mm X 25.4mm X 25.4mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5		

* Color does not apply to clear vented pattern