


# MILTEX® HI-LIGHT

## BECK ROOT ELEVATORS

### Beck Root Tip notch design

Elliptical sharp, thin, rounded tip for ease in penetrating interdental space



**Luxating step on tip**  
Prevents slipping



**DELB1**

### Permanent weld in handle

Preserves instrument integrity and longevity

### Finger platform on handle

Provides better grip on instrument

New Miltex® Beck Root Elevators are made from the highest quality German stainless steel and feature a permanent weld in the handle. The Beck Root Elevator's finger platform on the handle provides a better grip on the instrument.

Beck Root Elevators feature a sharp, thin rounded tip with luxating step. They are designed for easy penetration in the interdental space as well as 180° access to luxate in four directions, two mesial and two distal.



easy to penetrate interdental space



luxate in 4 directions

REF	Description
DELB1	Beck Elevator 2.8 mm Straight Tip
DELB2	Beck Elevator 3.5 mm Straight Tip
DELB3	Beck Elevator 4 mm Straight Tip
DELB4	Beck Elevator 3.5 mm Curved Right Tip
DELB5	Beck Elevator 3.5 mm Curved Left Tip



**DELB2**



**DELB3**



**DELB4**



**DELB5**

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

Illustrations and content provide general description only and may be subject to change.

Miltex, Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra LifeSciences Corporation. Copyright © 2009 Integra LifeSciences Corporation