

MILTEX HI-LIGHT

LUXATING ELEVATORS

Straight or Curved Blade

Range of styles and sizes for most applications.
Can be resharpener.

Unique Handle Shape

Creates less force on tooth, reducing chance of tooth fractures.



One-Piece Stainless Steel Construction

Tip will not separate from handle.
Provides strength and durability.

Miltex® Luxating Elevators are designed to make extractions as trauma free as possible by cutting the periodontal ligament while easing the tooth from the socket. The result is less tissue damage that promotes quicker healing and fewer complications. It is used by rocking and cutting rather than lifting and prying like a conventional elevator. The patterns include both straight and curved blades ranging from 2mm to 5mm. The sharp blade can be resharpener.



DELLUXS2
2mm Straight



DELLUXC2
2mm Curved



DELLUXS3
3mm Straight



DELLUXC3
3mm Curved



DELLUXS4
4mm Straight



DELLUXC4
4mm Curved



DELLUXS5
5mm Straight



DELLUXC5
5mm Curved

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

