


The LumbarTrak™ is the new Lumbar Retraction System for R&B Redmond™. The system is flexible enough to be used for many lumbar and thoracic spinal procedures from microdiscectomies to multilevel fusions. The entire system is made of titanium which is lighter and more radiolucent than stainless steel. The titanium blades provide enhanced visualization while still maintaining maximum exposure to the surgical site.

01. Color-coded blades (by length) allow for easy identification by the surgical staff to shorten the procedure time. A coordinated depth gauge is available to aid in blade length selection.
02. Blades come in lengths up to 130mm for larger patients and three widths for multiple procedures.
03. Quick and Long Connectors are available to attach multiple blades to the retractor. The Quick Connector can be attached with minimal effort. The Long Connector is used to adjust the lateral positioning of a blade.
04. A Retractor Handle is available for those situations where manual retraction is needed.

01






Retractors/Connectors:


	Description	Catalog#
01.	LumbarTrak™ Square Retractor	RB2601
02.	LumbarTrak™ Scoville Retractor	RB2602
03.	LumbarTrak™ Quick Connector	RB2603
04.	LumbarTrak™ Long Connector	RB2604





Blades:

Length	Color	Thin 24mm	Medium 45mm	Wide 60mm
30mm	Purple	RB2610	-----	-----
40mm	Light Blue	RB2611	-----	-----
50mm	Green	RB2612	RB2621	RB2620
60mm	Pink	RB2613	RB2623	RB2622
70mm	Dark Blue	RB2614	RB2625	RB2624
80mm	Gray	-----	RB2627	RB2626
90mm	Gold	RB2615	RB2629	RB2628
110mm	Purple	RB2616	RB2631	RB2630
130mm	Light Blue	RB2617	RB2633	RB2632





Right Blade 06.



Left Blade 07.



Right Blade 08.



Hooks/Sacral Blades:



Left Blade 09.


	Description	Catalog#	Specifications
01.	LumbarTrak™ Hook	RB2644	92mm
02.	LumbarTrak™ Hook	RB2643	72mm
03.	LumbarTrak™ Hook	RB2642	62mm
04.	LumbarTrak™ Hook	RB2641	52mm
05.	LumbarTrak™ Hook	RB2640	42mm
06.	LumbarTrak™ Sacral Blade Right	RB2647	40x30mm
07.	LumbarTrak™ Sacral Blade Left	RB2645	40x30mm
08.	LumbarTrak™ Sacral Blade Right	RB2648	45x60mm
09.	LumbarTrak™ Sacral Blade Left	RB2646	45x60mm



Handle/Gauge:




Description	Catalog#	Specifications
LumbarTrak™ Handle	RB2605	Blade not included
LumbarTrak™ Titanium Depth Gauge	RB2606	Color Coded
LumbarTrak™ Tray	RBT2600	




Complete Set: RB2600

Catalog#	Quantity	Catalog#	Quantity
RB2610	4	RB2628	2
RB2611	4	RB2629	2
RB2612	4	RB2630	2
RB2613	4	RB2631	2
RB2614	4	RB2632	2
RB2615	4	RB2633	2
RB2616	4	RB2640	2
RB2617	4	RB2641	2
RB2620	2	RB2642	2
RB2621	2	RB2643	2
RB2622	2	RB2644	2
RB2623	2	RB2645	2
RB2624	2	RB2646	2
RB2625	2	RB2647	2
RB2626	2	RB2648	2
RB2627	2	RB2601	1



Catalog#	Quantity
RB2602	1
RB2603	4
RB2604	4
RB2605	2
RB2606	1
RBT2600	1



Complete Set: RB2600



Surgical Instruments

USA and Canada
877-468-5572
609-750-4357 (Fax)

J. Jamner Surgical Instruments
9 Skyline Drive • Hawthorne, NY 10532

Outside USA and Canada
914-592-9050
914-592-8056 (Fax)

R&B Redmond and LumbarTrak are trademarks
of J. Jamner Surgical Instruments, Inc.
The Integra wave logo is a trademark of Integra LifeSciences Corporation.
© 2006 Integra LifeSciences Corporation. All Rights Reserved.



Printed in USA 10k NS1259-10/06