

Peapod IVD Rongeurs

Unique angled-up tip to reach into difficult positions

Cup sizes from 2mm to 4mm

Shaft lengths in 5 ½" and 7"




Ferris-Smith Rongeur

Ferris-Smith handle allows for greater power while maintaining control

Cup sizes from 2 x 10mm to 6 x 12mm

7" shaft length



Decker Forceps

Curved left and right for unique reaches

Cup size of 1.5 x 5mm

6" shaft length



Jackson IVD

70° up cup for extended reach

4 x 8mm cup size

7" shaft length



Williams Style Micro Disc IVD with 1 x 2 Teeth

1 x 2 teeth on end of cup for more secure grasping power

Shaft lengths in 5 ½" to 10"



IVD SIZING CHART



Angled-Up



Straight




Angled-Down

IVD Rongeurs Cup Size

1.5×10mm	2×10mm	2×12mm	2.5×7mm	3×7mm	3×10mm

3×12mm	4×7mm	4×8mm	4×10mm	4×14mm	5×10mm

5×14mm	6×10mm	6×12mm



R&BTM
REDMOND
Surgical Instruments

J. Jamner Surgical Instruments
9 Skyline Drive • Hawthorne, NY 10532

USA and Canada
 877-468-5572
 609-750-4357 (Fax)

Outside USA and Canada
 914-592-9050
 914-592-8056 (Fax)

R&B Redmond and Ultra-Coat are trademarks of J. Jamner Surgical Instruments, Inc.
 The Integra wave logo is a trademark of Integra LifeSciences Corporation.
 ©2006 J. Jamner Surgical Instruments, Inc. All Rights Reserved.

NS1270 - 09/06 Printed in USA 10k



REDMOND

Surgical Instruments

IVD RONGEURS



R&B Redmond™ offers a wide selection of inter-vertebral disc (IVD) rongeurs. From a bayonet style ideal for micro surgeries to ones with curved left and right jaws, R&B Redmond™ IVD Rongeurs can meet any need of the spinal surgeon.

Standard IVD Rongeurs

Available in cup sizes that range from 1.5 x 10mm to 5 x 10mm

Shaft lengths from 6" to 12"

Straight, angled-up and angled-down jaws allow access to multiple anatomic structures

Available with teeth in 7" sizes



Bayonet IVD Rongeurs


Bayonet design allows the user to work down a channel during microsurgery while keeping their hand out of the field of view

Available in Ultra-Coat™ coating, which is an extremely hard material to maintain sharpness for a longer period

7" shaft length, 12cm (4 3/4") working length

Available in cup sizes that range from 1.5 x 10mm to 4 x 10mm

Straight, angled-up and angled-down jaws



Caspar Type IVD Rongeur

Teeth at back of cup to allow for better grip

Cup sizes from 2 x 12mm to 6 x 16mm

Straight, angled-up and angled-down jaws

Shaft lengths in 6 1/4" and 7 1/4"



Wilde IVD Rongeur

Fenestrated cup to grab difficult tissue

Cup sizes from 2.5mm to 5mm

Straight, angled-up and angled-down jaws

Shaft lengths in 7" and 9"



Ordering Information

Standard IVD Rongeurs

	Length	Cushing 1.5 x 10mm	Love-Gruenwald 2 x 10mm	Spurling 3 x 10mm	Spurling 4 x 10mm	Spurling 5 x 10mm
Straight	6"		RB6201	RB6216	RB6231	
Up	6"		RB6202	RB6217	RB6232	
Down	6"		RB6203	RB6218	RB6233	
Straight	7"	RB6160	RB6204	RB6219	RB6234	
Up	7"	RB6161	RB6205	RB6220	RB6235	
Down	7"	RB6162	RB6206	RB6221	RB6236	
Straight	9"	RB6163	RB6207	RB6222	RB6237	RB6255
Up	9"	RB6164	RB6208	RB6223	RB6238	RB6256
Down	9"	RB6165	RB6209	RB6224	RB6239	RB6257
Straight	10"		RB6246	RB6249	RB6252	
Up	10"		RB6247	RB6250	RB6253	
Down	10"		RB6248	RB6251	RB6254	
Straight	11"		RB6210	RB6225	RB6240	RB6258
Up	11"		RB6211	RB6226	RB6241	RB6259
Down	11"		RB6212	RB6227	RB6242	RB6260
Straight	12"		RB6213	RB6228	RB6243	
Up	12"		RB6214	RB6229	RB6244	
Down	12"		RB6215	RB6230	RB6245	
Str w/ Teeth	7"		RB6301	RB6304	RB6307	
Up w/ Teeth	7"		RB6302	RB6305	RB6308	
Dwn w/ Teeth	7"		RB6303	RB6306	RB6309	

Stainless Bayonet IVD Rongeurs

	Working Length	Cushing 1.5 x 10mm	Cushing 2 x 10mm	Love-Gruenwald 3 x 10mm	Spurling 4 x 10mm
Straight	12cm	RB6601	RB6604	RB6621	RB6641
Up	12cm	RB6602	RB6605	RB6622	RB6642
Down	12cm	RB6603	RB6606	RB6623	RB6643

Ultra-Coat™ Bayonet IVD Rongeurs

	Working Length	Cushing 1.5 x 10mm	Cushing 2 x 10mm	Love-Gruenwald 3 x 10mm	Spurling 4 x 10mm
Straight	12cm	RB6801	RB6804	RB6821	RB6841
Up	12cm	RB6802	RB6805	RB6822	RB6842
Down	12cm	RB6803	RB6806	RB6823	RB6843

Caspar IVD Rongeurs

	Length	2 x 12mm	3 x 12mm	4 x 14mm	5 x 14mm	6 x 16mm
Straight	6 1/4"	RB6401	RB6402	RB6403	RB6404	RB6405
Up	6 1/4"		RB6406	RB6407	RB6408	
Down	6 1/4"		RB6409	RB6410		
Straight	7 1/4"	RB6411	RB6412	RB6413	RB6414	RB6415
Up	7 1/4"		RB6416	RB6417	RB6418	
Down	7 1/4"		RB6419	RB6420		

Ordering Information

Wilde IVD Rongeurs

	Length	2.5mm	3mm	4mm	5mm
Straight	7"	RB6007	RB6008	RB6009	RB6010
Up	7"	RB6011	RB6012	RB6013	RB6014
Straight	9"	RB6015	RB6016	RB6017	RB6018

Peapod IVD Rongeurs

	Length	2mm	2.5mm	3mm	4mm
Up	5 ½"	RB6126	RB6127	RB6128	
Up	7"		RB6129	RB6130	RB6131

Ferris Smith IVD Rongeurs

	Length	2 x 10mm	3 x 10mm	4 x 8mm	4 x 10mm	6 x 10mm	6 x 12mm
Straight	7"	RB6135	RB6136	RB6137	RB6138	RB6139	RB6140
Up	7"	RB6141	RB6142		RB6143	RB6144	
Down	7"	RB6145	RB6146		RB6147	RB6148	

Decker Forceps

	Length	1.5 x 5mm
Straight	6"	RB6005
Up	6"	RB6004
Down	6"	RB6003
Left	6"	RB6002
Right	6"	RB6001

Jackson IVD

	Length	4 x 8mm
70° Up	7"	RB6151

Williams Style Micro Disc IVD Rongeurs with 1 x 2 Teeth

	Length	1 x 2 Teeth
Straight	5 ½"	RB6152
Straight	7"	RB6154
Straight	10"	RB6155

