



LAPAROSCOPIC & BARIATRIC CASES

Make storage of endoscopic instruments easy!

PROTECT YOUR INSTRUMENTS

Now
available
in two sizes
to fit all of your
endoscopic
instruments




Protect your Laparoscopic Instruments with the JARIT® Laparoscopic Case and Laparoscopic Bariatric Case

Developed with the help of surgical technologies,
the case assures that instruments stay in place.

- The Laparoscopic Case fits JARIT® endoscopic instruments (32 cm and 37 cm).
- The Laparoscopic Bariatric Case fits JARIT® 45cm endoscopic instruments.

Bottom tray for storage of
ancillary instruments

Highly perforated
with rounded edges



Angled design promotes proper
drainage of lap instrument

Nylon coated
handle guide
helps to prevent
scratching
of insulated
instruments

PRODUCT ORDERING INFORMATION

The Laparoscopic Bariatric Case
Length 68.1 cm
Part # 600-620

The Laparoscopic Case
Length 57.9 cm
Part # 600-618



J. Jammer Surgical Instruments, Inc.
9 Skyline Drive, Hawthorne, New York 10532
www.jarit.com

USA and Canada
877-468-5572
609-750-4257 (Fax)

Outside USA and Canada
914-592-9050
914-592-8056 (Fax)

Jarrit is a registered trademark of J. Jammer Surgical Instruments, Inc., a company owned by Integra LifeSciences Corporation.
The Integra wave logo is a trademark of Integra LifeSciences Corporation.
© 2006 Integra LifeSciences Corporation. All Rights Reserved.

NS813 (12/05)
Printed in USA 10K