

Integra™

Neuro Spine Kerrison Rongeurs


Limits uncertainty with
durability and quality



INTEGRA™
LIMIT UNCERTAINTY

Surface coating options

All Ultra-Coat™ and Diamond™ Kerrison Rongeurs are sharpness warranted for a period of two (2) years under normal use. Should sharpening become necessary during the warranty period, Ultra-Coat and Diamond Kerrisons will be resharpened at no charge.



The Ultra-Coat kerrison rongeurs have a titanium nitride coating that more than quadruples the surface hardness vs. regular stainless steel, which extends the life of the cutting edge.



The Diamond kerrison rongeurs have a black ceramic coating that more than quadruples the surface hardness vs. regular stainless steel, which extends the life of the cutting edge and reduces intraoperative glare.



The Stealth kerrison rongeurs have a grey coating that improves the “glide” of the instrument. The Stealth™ coating does not increase the hardness of the cutting edge.

Detach Coated Rongeurs

The instrument disassembles to make it easier to clean and visually inspect.

510(k) cleared



Extended Length Kerrison Rongeurs

Additional lengths (to 13") addresses increasing patient BMI's. Available in all three surface coatings



Detach Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB586oDT	Ultra-Coat, 8", 40 degree up, 1mm
RB5862DT	Ultra-Coat, 8", 40 degree up, 2mm
RB5863DT	Ultra-Coat, 8", 40 degree up, 3mm
RB5864DT	Ultra-Coat, 8", 40 degree up, 4mm
RB5865DT	Ultra-Coat, 8", 40 degree up, 5mm
RD4891DT	Diamond, 8", 40 degree up, 1mm
RD4882DT	Diamond, 8", 40 degree up, 2mm
RD4853DT	Diamond, 8", 40 degree up, 3mm
RD4854DT	Diamond, 8", 40 degree up, 4mm
RD4855DT	Diamond, 8", 40 degree up, 5mm
RB4891DT	Stealth, 8", 40 degree up, 1mm
RB4882DT	Stealth, 8", 40 degree up, 2mm
RB4853DT	Stealth, 8", 40 degree up, 3mm
RB4854DT	Stealth, 8", 40 degree up, 4mm
RB4855DT	Stealth, 8", 40 degree up, 5mm

Detach Medium Kerrison Rongeurs - 15mm jaw opening

Catalog Number	Description
RB586oDTM	Ultra-Coat, 8", 40 degree up, 1mm
RB5862DTM	Ultra-Coat, 8", 40 degree up, 2mm
RB5863DTM	Ultra-Coat, 8", 40 degree up, 3mm
RB5864DTM	Ultra-Coat, 8", 40 degree up, 4mm
RB5865DTM	Ultra-Coat, 8", 40 degree up, 5mm
RD4821DTM	Diamond, 8", 40 degree up, 1mm
RD4822DTM	Diamond, 8", 40 degree up, 2mm
RD4823DTM	Diamond, 8", 40 degree up, 3mm
RD4824DTM	Diamond, 8", 40 degree up, 4mm
RD4825DTM	Diamond, 8", 40 degree up, 5mm
RB4821DTM	Stealth, 8", 40 degree up, 1mm
RB4822DTM	Stealth, 8", 40 degree up, 2mm
RB4823DTM	Stealth, 8", 40 degree up, 3mm
RB4824DTM	Stealth, 8", 40 degree up, 4mm
RB4825DTM	Stealth, 8", 40 degree up, 5mm

Extended Length Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB5362	Ultra-Coat, 13", 40 degree up, 2mm
RB5363	Ultra-Coat, 13", 40 degree up, 3mm
RB5364	Ultra-Coat, 13", 40 degree up, 4mm
RB5365	Ultra-Coat, 13", 40 degree up, 5mm
RB5366	Ultra-Coat, 13", 40 degree up, 6mm
RD41192	Diamond, 11", 40 degree up, 2mm
RD41193	Diamond 11", 40 degree up, 3mm
RD41194	Diamond 11", 40 degree up, 4mm
RD41195	Diamond 11", 40 degree up, 5mm
RD41196	Diamond 11", 40 degree up, 6mm

Medium Kerrison Rongeurs - 15mm jaw opening

Catalog Number	Description
RB4721M	Stealth, 7", 40 degree up, 1mm
RB4751M	Stealth, 7", 40 degree up, 1.5mm
RB4722M	Stealth, 7", 40 degree up, 2mm
RB4723M	Stealth, 7", 40 degree up, 3mm
RB4724M	Stealth, 7", 40 degree up, 4mm
RB4725M	Stealth, 7", 40 degree up, 5mm
RB4726M	Stealth, 7", 40 degree up, 6mm
RB4821M	Stealth, 8", 40 degree up, 1mm
RB4851M	Stealth, 8", 40 degree up, 1.5mm
RB4822M	Stealth, 8", 40 degree up, 2mm
RB4823M	Stealth, 8", 40 degree up, 3mm
RB4824M	Stealth, 8", 40 degree up, 4mm
RB4825M	Stealth, 8", 40 degree up, 5mm
RB4826M	Stealth, 8", 40 degree up, 6mm
RD4821M	Diamond, 8", 40 degree up, 1mm
RD4851M	Diamond, 8", 40 degree up, 1.5mm
RD4822M	Diamond, 8", 40 degree up, 2mm
RD4823M	Diamond, 8", 40 degree up, 3mm
RD4824M	Diamond, 8", 40 degree up, 4mm
RD4825M	Diamond, 8", 40 degree up, 5mm
RD4826M	Diamond, 8", 40 degree up, 6mm
RB586oM	Ultra-Coat, 8", 40 degree up, 1mm
RB5862M	Ultra-Coat, 8", 40 degree up, 2mm
RB5863M	Ultra-Coat, 8", 40 degree up, 3mm
RB5864M	Ultra-Coat, 8", 40 degree up, 4mm
RB5865M	Ultra-Coat, 8", 40 degree up, 5mm

Extended Length Bayonet Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB5662	Ultra-Coat, Bayonet 11", 40 degree up, 2mm
RB5663	Ultra-Coat, Bayonet 11", 40 degree up, 3mm
RB5664	Ultra-Coat, Bayonet 11", 40 degree up, 4mm
RB5665	Ultra-Coat, Bayonet 11", 40 degree up, 5mm
RB5666	Ultra-Coat, Bayonet 11", 40 degree up, 6mm

For a complete listing of kerrison rongeurs refer to the catalog

For more information or to place an order, please contact:
 Integra ■ 311 Enterprise Drive, Plainsboro, NJ 08536
 USA and Canada: 877.468.5572 ■ 914.592.9050 (Outside USA)
 800.275.7192 (Fax)
integralife.com

