

Miltex
Redefining Excellence



STERILIZATION CONTAINERS
AND SPD SUPPLIES



INTRODUCTION

Miltex has solicited the help and experience of our Customers, Surgical Dealer Partners and the finest European Manufacturers to design and construct one of the best Sterilization Containers available today.

We have included all of the necessary choices while streamlining the selection to offer not only an excellent Container but one that is economical and well suited for your facility budget and instrument sterilization and storage requirements.

Miltex® – Redefining Excellence


MILTEX® STERILIZATION CONTAINER SYSTEM

Miltex® Containers.....Designed and Constructed for the Customer

- Container lid and bottom are both of one-piece construction and have a uniform, anodized aluminum alloy finish
- Lightweight, yet durable enough to accommodate acceptable weight loads of surgical instrumentation
- Handles stop at 90° to aid lifting; low profile latch allows simple and secure lid closing
- Retention plate employs an easy "click and lock" closure over disposable, quick-to-change filter
- Dual silicone gaskets within retention plate and around perimeter ensure complete Container and filter seal
- Basket offering includes both solid wire mesh and perforated metal designs to allow for maximum ultrasonic cleaning; handles can be positioned to assist with secure and comfortable grasp for removal under sterile conditions



SUMMARY OF MILTEX® STERILIZATION CONTAINER PRODUCT FEATURES



ONE-HALF SIZE (11-1/4" X 11") CONTAINER WITH LID

- Vacuum Sterilizer Containers - Solid Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5100-10	3-5110-10	4"	1/2-Size Complete Container, Silver
3-5100-13	3-5110-13	5-1/3"	1/2-Size Complete Container, Silver
3-5100-15	3-5110-15	6"	1/2-Size Complete Container, Silver

*Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.*



COLORED LIDS (bottoms sold separately)

**Includes Filter Retention
Plate and Latches**

REF	Description
3-6100-00	1/2-Size Container Lid Only, Silver
3-6100-01	1/2-Size Container Lid Only, Blue
3-6100-02	1/2-Size Container Lid Only, Red
3-6100-03	1/2-Size Container Lid Only, Gold
3-6100-04	1/2-Size Container Lid Only, Green
3-6100-05	1/2-Size Container Lid Only, Black



CONTAINER BOTTOMS VACUUM/GRAVITY (lids sold separately)

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5200-10	3-5210-10	4"	1/2-Size Container Bottom Only
3-5200-13	3-5210-13	5-1/3"	1/2-Size Container Bottom Only
3-5200-15	3-5210-15	6"	1/2-Size Container Bottom Only

3-5210-10



3-5200-13



3-5200-15



ONE-HALF SIZE WIRE/PERFORATED BASKETS

**Wire Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5172-05	1/2-Size Small (2" approx.)	245 x 255 x 50
3-5172-07	1/2-Size Medium (3" approx.)	245 x 255 x 70
3-5172-10	1/2-Size Large (4" approx.)	245 x 255 x 100

**Perforated Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5174-05	1/2-Size Small (2" approx.)	240 x 240 x 50
3-5174-07	1/2-Size Medium (3" approx.)	240 x 240 x 70
3-5174-09	1/2-Size Large (4" approx.)	240 x 240 x 90

**Perforated Basket with Feet**

Feet help to raise basket over retention plate and are ideal for use with gravity perforated bottom containers. Feet are removable.

REF	Description	Dimensions in mm (LxWxH)
3-5174-05WF	1/2-Size Small (2" approx.)	242 x 252 x 59
3-5174-07WF	1/2-Size Medium (3" approx.)	242 x 252 x 78
3-5174-09WF	1/2-Size Large (4" approx.)	242 x 252 x 108

THREE-QUARTER SIZE (18-1/4" X 11") CONTAINER WITH LID

- Vacuum Sterilizer Containers - Solid Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5300-10	3-5310-10	4"	3/4-Size Complete Container, Silver
3-5300-13	3-5310-13	5-1/3"	3/4-Size Complete Container, Silver
3-5300-15	3-5310-15	6"	3/4-Size Complete Container, Silver

*Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.*



COLORED LIDS (bottoms sold separately)

Includes Filter Retention Plate and Latches

REF	Description
3-6300-00	3/4-Size Container Lid Only, Silver
3-6300-01	3/4-Size Container Lid Only, Blue
3-6300-02	3/4-Size Container Lid Only, Red
3-6300-03	3/4-Size Container Lid Only, Gold
3-6300-04	3/4-Size Container Lid Only, Green
3-6300-05	3/4-Size Container Lid Only, Black



CONTAINER BOTTOMS VACUUM/GRAVITY (lids sold separately)

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5400-10	3-5410-10	4"	3/4-Size Container Bottom Only
3-5400-13	3-5410-13	5-1/3"	3/4-Size Container Bottom Only
3-5400-15	3-5410-15	6"	3/4-Size Container Bottom Only

3-5400-10



3-5410-13



3-5400-15



THREE-QUARTER SIZE WIRE/PERFORATED BASKETS

**Wire Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5372-05	3/4-Size Small (2" approx.)	405 x 255 x 50
3-5372-07	3/4-Size Medium (3" approx.)	405 x 255 x 70
3-5372-10	3/4-Size Large (4" approx.)	405 x 255 x 100

**Perforated Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5374-05	3/4-Size Small (2" approx.)	400 x 240 x 50
3-5374-07	3/4-Size Medium (3" approx.)	400 x 240 x 70
3-5374-09	3/4-Size Large (4" approx.)	400 x 240 x 90

**Perforated Basket with Feet**

Feet help to raise basket over retention plate and are ideal for use with gravity perforated bottom containers. Feet are removable.

REF	Description	Dimensions in mm (LxWxH)
3-5374-05WF	3/4-Size Small (2" approx.)	402 x 252 x 59
3-5374-07WF	3/4-Size Medium (3" approx.)	402 x 252 x 78
3-5374-09WF	3/4-Size Large (4" approx.)	402 x 252 x 108

FULL-SIZE (23" X 11") CONTAINER WITH LID

- Vacuum Sterilizer Containers - Solid Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5500-10	3-5510-10	4"	Full-Size Complete Container, Silver
3-5500-13	3-5510-13	5-1/3"	Full-Size Complete Container, Silver
3-5500-15	3-5510-15	6"	Full-Size Complete Container, Silver

Note: Complete container assembly only available in silver.

See below for other color options with lids and bottoms sold separately.



COLORED LIDS (bottoms sold separately)

Includes Filter Retention Plate and Latches

REF	Description
3-6500-00	Full-Size Container Lid only, Silver
3-6500-01	Full-Size Container Lid only, Blue
3-6500-02	Full-Size Container Lid only, Red
3-6500-03	Full-Size Container Lid only, Gold
3-6500-04	Full-Size Container Lid only, Green
3-6500-05	Full-Size Container Lid only, Black



CONTAINER BOTTOMS VACUUM/GRAVITY (lids sold separately)

Vacuum Sterilizer Containers **Gravity Sterilizer Containers**

REF	REF	Height	Description
3-5600-10	3-5610-10	4"	Full-Size Container Bottom Only
3-5600-13	3-5610-13	5-1/3"	Full-Size Container Bottom Only
3-5600-15	3-5610-15	6"	Full-Size Container Bottom Only



FULL-SIZE WIRE/PERFORATED BASKETS

**Wire Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5572-05	Full-Size Small (2" approx.)	540 x 255 x 50
3-5572-07	Full-Size Medium (3" approx.)	540 x 255 x 70
3-5572-10	Full-Size Large (4" approx.)	540 x 255 x 100

**Perforated Basket**

REF	Description	Dimensions in mm (LxWxH)
3-5574-05	Full-Size Small (2" approx.)	540 x 240 x 50
3-5574-07	Full-Size Medium (3" approx.)	540 x 240 x 70
3-5574-09	Full-Size Large (4" approx.)	540 x 240 x 90

**Perforated Basket with Feet**

Feet help to raise basket over retention plate and are ideal for use with gravity perforated bottom containers. Feet are removable.

REF	Description	Dimensions in mm (LxWxH)
3-5574-05WF	Full-Size Small (2" approx.)	538 x 252 x 59
3-5574-07WF	Full-Size Medium (3" approx.)	538 x 252 x 78
3-5574-09WF	Full-Size Large (4" approx.)	538 x 252 x 108

TAMPER EVIDENT LOCKS



REF	Description
3-5935-00	Blue, box of 1,000

IDENTIFICATION TAGS

REF	Description
3-5940-01	ID Tag, silver
3-5940-02	ID Tag, blue
3-5940-03	ID Tag, red
3-5940-04	ID Tag, gold
3-5940-05	ID Tag, green
3-5940-06	ID Tag, black



Contents/Delivery Address - Custom engraving available with order.
Each tag allows up to 14 characters.

INDICATOR CARDS



REF	Description
3-5930-00	Indicates for use in gas or steam, box of 250

FILTERS



REF	Description
3-5900-15	Single use, round 7-1/2", box of 1,000

MILTEX® CARE & CLEANING PRODUCTS

See page 16-17 for more product details



ENDOSCOPIC INSTRUMENT RACK

For The Miltex® Sterilization Container System



REF	Description
3-5770-00	Endoscopic Instrument Rack 19-3/4" L (50.2 cm) x 8-3/4" W (22.2 cm) x 4" H (10.2 cm)



FEATURES AND BENEFITS

- Compatible with 6" Full-Size Sterilization Container.
- Rack fits inside Full-Size perforated or wire basket and allows ample storage underneath for other instruments and auxiliary items.
- Holds up to 10 instruments and suitable for Arthroscopic, Laparoscopic and Endoscopic instruments ranging from 3mm to 8mm OD.
- Latex-Free and constructed with medical grade silicone grips and stainless steel frame.
- Cushioned grips protect delicate insulation of endoscopic instruments and provides protection during transport.
- Offset side hand grips facilitate easy and aseptic removal of rack and instruments.



3-70 to 3-74



3-75 to 3-80



3-81 to 3-85



3-86 to 3-89



3-90 to 3-95



3-96 to 3-99

INSTRUMENT STRINGERS



REF	Description	Size
3-70	Stringer - Side Lock, Wide	2-1/2" x 8"
3-71	Stringer - Side Lock, Wide	2-1/2" x 10"
3-72	Stringer - Side Lock, Wide	2-1/2" x 12"
3-74	Stringer - Side Lock, Wide	2-1/2" x 16"
3-75	Stringer - Side Lock	2-1/2" x 6"
3-76	Stringer - Side Lock	2-1/2" x 8"
3-77	Stringer - Side Lock	2-1/2" x 10"
3-78	Stringer - Side Lock	2-1/2" x 12"
3-79	Stringer - Side Lock	2-1/2" x 16"
3-80	Stringer - Side Lock	2-1/2" x 18"
3-81	Stringer - Side Lock	5" x 8"
3-82	Stringer - Side Lock	5" x 10"
3-83	Stringer - Side Lock	5" x 12"
3-84	Stringer - Side Lock	5" x 16"
3-85	Stringer - Side Lock	5" x 18"
3-86	Stringer - Open U	2-1/2" x 6"
3-87	Stringer - Open U	2-1/2" x 8"
3-88	Stringer - Open U	2-1/2" x 10"
3-89	Stringer - Open U	2-1/2" x 12"
3-90	Stringer - Center Lock	2-1/2" x 6"
3-91	Stringer - Center Lock	2-1/2" x 8"
3-92	Stringer - Center Lock	2-1/2" x 10"
3-93	Stringer - Center Lock	2-1/2" x 12"
3-94	Stringer - Center Lock	2-1/2" x 16"
3-95	Stringer - Center Lock	2-1/2" x 18"
3-96	Stringer - Closed U	2-1/2" x 6"
3-97	Stringer - Closed U	2-1/2" x 8"
3-98	Stringer - Closed U	2-1/2" x 10"
3-99	Stringer - Closed U	2-1/2" x 12"

Miltex® Instrument Storage System featuring the “Tuck-Away™” Lid



1. Unique flow-through design.
2. Medical grade silicone strip provides instrument retention when lid is closed.
3. Full length hinge provides stable and secure attachment to cassette.
4. Quick and easy push button release provides a positive lock when closed.
5. Pre-stamped tabs allow for custom coding of cassettes.
6. Double-hinged lid folds under cassette for easy storage during use.
7. Self-lubricating plastic hinge pin will not corrode and provides smooth mechanical movement.
8. Passivated stainless steel provides corrosion protection and an aesthetic finish.

UNIVERSAL CASSETTE WITHOUT MAT


REF	Size
4-6835	6-1/2 x 6-1/4 x 1-3/16
4-083000	8 x 3-13/16 x 1-3/16
4-084000	8 x 4-7/8 x 1-3/16
4-085000	8 x 6 x 1-3/16
4-089100	8 x 9 x 1-3/16
4-008122	8 x 11 x 1-3/16
4-008115	8 x 11 x 2
4-009126	9-1/8 x 13 x 1-3/16
4-6815	15 x 8-1/4 x 1-3/16

SLIMLINE SERIES WITHOUT MAT

3-072007	7 x 2-5/8 x 11/16*
3-082008	8 x 2-5/8 x 11/16*
STD5L27	7 x 2-5/8 x 1*
3-072014	7 x 4-5/8 x 11/16*

*Single-hinge design

Latex-Free Silicone Mats



1. The raised mat (or feet) assures adhesion of the mat to the tray surface.
2. Holes positioned on mat for proper drainage.
3. Soft cushion fingers protect your instrument investment.
4. Customized to fit perfectly into Miltex® Cassettes.

MAT ONLY

REF	Size	Description
4-6835MAT	6-1/2 x 6-1/4	for Cassette 4-6835
4-083000MAT	8 x 3-13/16	for Cassette 4-083000
4-084000MAT	8 x 4-7/8	for Cassette 4-084000
4-085000MAT	8 x 6	for Cassette 4-085000
4-089100MAT	8 x 9	for Cassette 4-089100
4-008122MAT	8 x 11	for Cassette 4-008122 & 4-008115
4-009126MAT	9-1/8 x 13	for Cassette 4-009126
4-6815MAT	15 x 8-1/4	for Cassette 4-6815

MILTEX® Tray with Silicone Mat


3-200100 7" x 11" (17.8 x 27.9 cm), single layer



3-200100

MILTEX® Tray with Silicone Mats

3-200200 7" x 11" (17.8 x 27.9 cm), double layer



3-200200

MILTEX® Tray with Silicone Mat

3-200300 4" x 8" (10.2 x 20.3 cm), single layer



3-200300

Small Plastic Tray with Silicone Mat

PM-9323* 2-3/4" x 1-3/4" x 3/4" (73 x 44 x 19 mm)



PM-9323

Microsurgical Plastic Tray with Silicone Mat

PM-9326* 10" x 5-3/4" x 3/4" (254 x 146 x 19 mm)



PM-9326

Microsurgical Plastic Tray with Silicone Mat

PM-9325* 11" x 7" x 3/4" (279 x 178 x 19 mm)



PM-9325

MILTEX® INSTRUMENT CARE & CLEANING PRODUCTS

The MILTEX® care and cleaning product portfolio is designed to provide optimum cleaning performance and will help to prolong the service life of valuable surgical instrumentation.

MILTEX® Instrument Prep Enzyme Foam

Ready-to-use foaming spray for pre-cleaning of soiled instruments and scopes (rigid or flexible). Ideal for applications where instrument cleaning is not immediately available. Simply spray on instruments until ready for decontamination and cleaning.

3-760 24 oz. (0.71 liter) spray bottles (case of 12)*



3-760

MILTEX® Spray Lube (Surgical Instrument Lubricant)

Pre-mixed, non-silicone and anti-corrosive formula in a ready to use 8 oz. spray pump bottle. Steam permeable and water soluble. Helps prevent spotting, staining and rusting. Lubricates instrument box locks, screw locks, scissors blades and any other metal on metal moving parts.

3-700 8 fl. oz. (0.24 liter) bottles (case of 12)*



3-700

MILTEX® EZ-ZYME® All Purpose Enzyme Cleaner

Ideal for cleaning hard to remove proteinaceous (mucus & tissue) and lipid (fatty) materials from the surfaces of surgical and dental instruments, linen, water submersible equipment, or where mechanical scrubbing is impractical.

- Easy to Use
- Neutral pH
- Multi-use as an instrument presoak, ultrasonic cleaner, laundry presoak and evacuation system cleaner
- Space-saving or bulk packaging options

3-750 3/4 oz. (22 ml) packets (dispenser box of 32)
Contents of 1-box makes 32 gallons of solution

3-755 1-gallon bottle (3.79 liter)
Contents of 1-gallon container makes 170 gallons of solution



3-755



3-750

MILTEX® Surgical Instrument Cleaner

Super concentrate formula for manual and ultrasonic cleaning or soaking of surgical instruments, accessories and glassware.

- Neutral pH
- Biodegradable
- Rinses free, leaves no residue
- Phosphate-free
- Concentrated

3-720 8 fl. oz. (0.24 liter) bottles (case of 12)*
Contents of one 8 oz. bottle makes 32 gallons of solution

3-725 1-gallon bottle (3.79 liter)
Contents of 1-gallon container makes 512 gallons of solution

3-726 1/4 fl. oz. (7.4 ml) packets (dispenser box of 80)
Contents of 1-box makes 80 gallons of solution



3-725



3-720



3-726

*(Contact Distributor for availability of individual units)

MILTEX® Surgical Instrument Cleaner & Lubricant

A dual-purpose instrument cleaner and lubricant. Chemically engineered for use as a presoak or for manual cleaning, washer/decontaminators, ultrasonic cleaners and cart washers.

- Cleaner & Brighter Instruments
- Removes Stains, Blood and Mineral Deposits
- Biodegradable
- Concentrate, Lower Cost Per Use

3-710 1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

**3-710****MILTEX® Instrument Lubricant**

A water emulsible soaking agent works to extend the life of surgical instruments. Concentrated non-silicone lubricant and rust inhibitor. Use daily to keep moving parts operating smoothly.

- Formulated to extend the working life of surgical instruments
- Helps Prevent Corrosion
- Formulated with gentle, non-aggressive ingredients
- Frees Up Box Locks - Lubricates moving parts of surgical instruments

3-705 1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

**3-705****MILTEX® Surgical Instrument Stain Remover**

Powder formula removes stains, tarnish and restores luster to surgical instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, etc. Ideal for use on stainless steel, chrome, nickel, brass, copper and silver.

- Effective on instruments, autoclaves, ultrasonic cleaners or any stainless steel, tile or plated surface
- Will not harm instrument surfaces
- Easy, fast, convenient to use
- Economical - a little goes a long way

3-740 3 oz. (85 grams), (case of 12)*

**3-740****MILTEX® Instrument Care Kit System**

A great way to try various MILTEX® products!

3-800 Spray Lube (3-700), Surgical Instrument Cleaner (3-720), Surgical Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon Brush and Stainless Steel Brush

**3-800**

*(Contact Distributor for availability of individual units)

MILTEX® Instrument Care Training DVD

25-minute Instructional DVD assists practitioners and staff to properly care for instruments and to keep them in top working condition. DVD offers step-by-step guidance on handling, cleaning, rinsing, inspecting, lubricating and sterilizing surgical instruments.

3-810 MILTEX® Instrument Care DVD



3-810

MILTEX® Nylon & Stainless Steel Brushes

- Nylon bristles for effective cleaning of instrument box locks
- Stainless steel bristles for cleaning serrations
- Autoclavable

3-1000 7-1/4" (18.4 cm), Nylon Bristles (Set of 3)

3-1001 7-1/4" (18.4 cm), Stainless Steel Bristles (Set of 3)

21-379 6-1/4" (15.9 cm), Stainless Wire Brush



SILITONE™ Color Bands

Silitone™ color bands are made from medical-grade silicone and can withstand temperatures up to 450° F. Silitone™ bands are expandable from 4 mm to 12 mm and are available in a variety of colors. 50 pieces per pack except kits.

REF	Description
VSTSA	Grey - 50 pieces
VSTSB	Blue - 50 pieces
VSTSD	Gold - 50 pieces
VSTSG	Green - 50 pieces
VSTSK	Black - 50 pieces
VSTSL	Light Blue - 50 pieces
VSTSN	Brown - 50 pieces
VSTSP	Pink - 50 pieces
VSTSR	Red - 50 pieces
VSTSU	Purple - 50 pieces
VSTSW	White - 50 pieces
VTSY	Yellow - 50 pieces
VST	Kit of 12 colors - 480 pieces (40 of each color)
VSTA	Kit of 12 colors - 120 pieces (10 of each color)



Will work with ring-handled instruments!



VST

(Exercise caution when applying color bands to instruments with sharp points.)

TAPE N' TELL™ Adhesive Strips

The Tape N' Tell™ adhesive strips are ideal for color coding instruments and procedural kits.

REF	Description
017-48220	TAPE N' TELL WHITE
017-48221	TAPE N' TELL YELLOW
017-48222	TAPE N' TELL RED
017-48223	TAPE N' TELL BLUE
017-48224	TAPE N' TELL GREEN
017-48225	TAPE N' TELL BLACK
017-48226	TAPE N' TELL ORANGE
017-48227	TAPE N' TELL BROWN
017-48230	TAPE N' TELL 8 ASSORTED



TIP-IT™ Instrument Guards

Tip-It™ Instrument Guards protect instrument tips and edges from damage during handling and sterilization.

- PVC medical grade plastic
- Latex-free material
- Large holes (vented styles) to reduce vacuum and moisture build up for easy removal from instrument
- Steam and gas sterilizable
- 50 per pack



REF Clear Vented	REF Tinted Vented	REF Solid Vented	REF Solid Non-Vented	Dimension	Size Code	Color
3-2501C	3-2501T	3-2501V	3-2501	1/16" x 3/4" (1.6 mm x 19 mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" x 3/4" (2 mm x 19 mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" x 3/4" (2.8 mm x 19 mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" x 1" (3.2 mm x 25.4 mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" x 1" (4.8 mm x 25.4 mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" x 3/4" (10.3 mm x 19 mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" x 3/8" x 1" (1.6 mm x 9.5 mm x 25.4 mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" x 5/8" x 1" (2 mm x 16 mm x 25.4 mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" x 1" x 1" (3.2 mm x 25.4 mm x 25.4 mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5	-	-

* Color does not apply to clear vented pattern

Sterilization Products

Miltex offers a wide range of sterilization products which can be used in conjunction with the Universal cassettes.

CSR WRAP


A suitable CSR wrap for steam, ETO gas and radiation sterilization processes. CSR wraps are electronically inspected for defects such as pin holes and are monitored for thickness consistency. CSR wraps are one of the lowest linting products available on the market.

Features include:

- Resists memory
- Allows faster steam penetration than muslin
- Repels water and OR fluids
- Low linting
- Meets established standards for tensile strength and tear resistance
- Excellent barrier to airborne bacteria
- Accepts sterilization tape and labels
- Biodegradable

Available in three sizes—15", 20", and 24" square

Item#	Dimensions	Case quantity
STD1515	15" x 15"	1000/case
STD2020	20" x 20"	500/case
STD2424	24" x 24"	500/case



1. CSR Wrap
2. Indicator Tapes

INDICATOR TAPES

Used to visually identify processed from unprocessed packs. Miltex® Indicator Tapes are manufactured from a specially treated crêpe paper and are coated with a high performance, pressure-sensitive adhesive that will not come undone when subjected to wet or dry heat.

Features include:

- *Adheres to any kind of material*
- *Distinctive color-change*
- *Accommodates pack expansion during sterilization*
- *Accepts labels and writing*

Supplied in 60 yard rolls, Miltex® Indicator Tapes are available in steam, vapor and dry heat styles.

Item#	Description
STDSTEAM	3/4" Steam, 60 yd.
STDCHEM	1/2" Chemical Vapor, 60 yd
STDDRY	1" Dry Heat, 60 yd.

STEAM INDICATOR STRIPS

Per CDC guidelines, a chemical indicator should be used with each load in conjunction with routine spore testing. The Miltex® Steam Indicator Strip turns brown when sterilization conditions have been met and is certified for use with all steam sterilization processes.

Packaged 250 strips per box, the Steam Indicator strips may be maintained as a permanent record of sterilization

Item#	Description
STDIND	Steam Indicator Strips, 250/bx

Surgical Instruments, Floor Grade Quality, Sterile

Single-Use Instruments

- Convenience and cost savings provide an exceptional value
- Eliminate the time and expense of re-processing soiled instruments
- Includes all of the most common instrument patterns for routine surgical needs
- Sterile, 2-year shelf life



REF	Description	
ST5-120	MAYO, 5-1/2", ST, B/B	50/Bx
ST5-124	MAYO, 6-3/4", ST, B/B	50/Bx
ST5-16	OR Scissors, 5-1/2", ST, S/B	50/Bx
ST5-180	METZENBAUM, 5-1/2", Curved	50/Bx
ST5-182	METZENBAUM, 7", Curved	50/Bx
ST5-304	IRIS Scissors, 4-1/2", Straight	50/Bx
ST5-306	IRIS Scissors, 4-1/2", Curved	50/Bx
ST5-504SS	LISTER Bandage Scissors, 5-1/2"	50/Bx
ST9-102	LITTAUER Stitch Scissors, 4-1/2"	50/Bx

SCISSORS

ST6-118	ADSON Forceps, 4-3/4"	50/Bx
ST6-120	ADSON Forceps, 4-3/4", 1X2	50/Bx
ST6-44	TISSUE Forceps, 5-1/2", 1X2	50/Bx
ST6-8	DRESSING Forceps, 5-1/2"	50/Bx
ST7-138	ROCHESTER-PEAN, 6-1/4", Curved	50/Bx
ST7-2	HALSTEAD Forceps, 5", Straight	50/Bx
ST7-4	HALSTEAD Forceps, 5", Curved	50/Bx
ST7-36	KELLY Forceps, 5-1/2", Straight	50/Bx
ST7-38	KELLY Forceps, 5-1/2", Curved	50/Bx
ST7-42	CRILE Forceps, 5-1/2", Straight	50/Bx
ST7-44	CRILE Forceps, 5-1/2", Curved	50/Bx

FORCEPS & HEMOSTATS

ST8-44	MAYO-HEGAR Needle Holder, 6"	50/Bx
ST8-50	CRILE-WOOD Needle Holder, 6"	50/Bx
ST8-6	WEBSTER Needle Holder, 5"	50/Bx

NEEDLE HOLDERS



SPONGE BOWLS

REF	Size
3-900	4-5/16" x 2-13/64" (12 oz)
3-901	4-1/4" x 2" (16 oz)
3-902	5-5/8" x 2-9/32" (22 oz)
3-903	6-27/64" x 2-1/4" (30 oz)
3-904	6-27/64" x 2-3/4" (36 oz)
3-905	7-9/16" x 3-3/16" (47 oz)
3-906	8-3/16" x 3-5/16" (64 oz)



IODINE & MEDICINE CUPS

REF	Description	Size
3-912	Iodine Cup	3-5/16" x 2" (6 oz)
3-913	Iodine Cup	5-1/32" x 2-23/64" (14 oz)
3-914	Medicine Cup, Graduated	2-5/32" x 2-5/32" (2 oz)



TUMBLERS

REF	Size
3-907	2-7/8 x 3-1/8 (7 oz)



SOLUTION BASINS

REF	Size
3-908	12-33/64" x 3-35/64" (4-1/2 qt)
3-909	13-25/32" x 4-17/32" (6 qt)
3-910	13-25/32" x 4-59/64" (7 qt)

Items in **bold** represent the most popular patterns.

**FORCEPS JARS**

REF	Size
3-915	1-19/64" x 3-5/32"
3-916	2-3/16" x 4-21/64"
3-917	2-3/16" x 7-5/64"

**DRESSING JARS WITH COVERS**

Cover included with jar

REF	Size
3-953	4" x 2-5/8" (1/2 qt)
3-954	5" x 6-3/8" (2 qt)
3-955	5-7/8" x 7-5/64" (3 qt)

**EMESIS BASINS**

REF	Size
3-918	7-7/8" x 3-3/4" x 1-19/64" (10 oz)
3-919	8-5/32" x 3-55/64" x 1-17/32" (12 oz)
3-920	9-23/32" x 4-51/64" x 1-11/16" (20 oz)

**INSTRUMENT TRAYS**

REF	Size
3-934	10" x 6-1/2" x 23/32", rolled edge
3-935	14" x 10" x 5/8", rolled edge
3-936	21-11/32" x 16-7/16" x 7/8", rolled edge

**PERFORATED MAYO TRAYS**

REF	Size
3-921	10" x 6-1/2" x 23/32"
3-922	14" x 10" x 5/8"
3-924	17" x 11-5/8" x 3/4"
3-925	19-11/64" x 12-23/32" x 33/64"

**NON-PERFORATED MAYO TRAYS**

REF	Size
3-926	10" x 6-1/2" x 23/32"
3-927	14" x 10" x 5/8"
3-929	17" x 11-5/8" x 3/4"
3-930	19-11/64" x 12-23/32" x 33/64"

Items in **bold** represent the most popular patterns.



INSTRUMENT TRAYS & COVERS cover sold separately

REF	Description	Size
3-945	Instrument Tray, 12, 2-1/2, Solid	12-23/32" x 10-27/64" x 2-1/2"
3-946	Instrument Tray, 12, 3-3/4, Solid	12-23/32" x 10-27/64" x 3-3/4"
3-947	Flat Cover, 12, For REF 3-945 & 3-946	12-29/32" x 10-35/64" x 7/16"
3-948	Instrument Tray, 10, 1, Solid	10-25/64" x 6-49/64" x 1-27/32"
3-949	Flat Cover, 10, For REF 3-948	10-15/32" x 6-61/64" x 7/16"



CATH TRAYS & COVERS cover sold separately

REF	Description	Size
3-937	Cath Tray, 8 x 5 x 1	8-25/32" x 4-61/64" x 1-49/64"
3-938	Cath Tray, 8-7/8 x 5 x 3	8-7/8" x 5" x 3"
3-939	Flat Cath Cover, 8, For REF 3-937 & 3-938	8" x 5" x 7/16"
3-940	Dome Cath Cover, 8, For REF 3-937 & 3-938	8" x 5" x 7/16"
3-941	Cath Tray, 8 x 3 x 1-1/2, Narrow	8-35/64" x 3-5/64" x 1-1/2"
3-942	Dome Cath Cover, 8, Narrow, For REF 3-941	8-21/32" x 3-5/16" x 19/32"
3-943	Cath Tray, 15-3/4 x 2-3/4 x 2, Long	15-3/4" x 2-3/4" x 2"
3-944	Dome Cath Cover, 15, Long, For REF 3-943	15-13/16" x 3-1/32" x 29/32"

Items in **bold** represent the most popular patterns.



LIMITED WARRANTY

Miltex hereby warrants all Miltex® Containers purchased directly from Miltex or an Authorized Dealer in accordance with the following terms and conditions:

The Miltex Container is guaranteed to be free from defects in material and workmanship when used for its intended purpose. Any Miltex Container that proves defective in workmanship or material will either be repaired or replaced, at Miltex discretion, without charge.

Containers not purchased from an authorized source may not be authentic or meet Miltex quality standards.

Container Gaskets Warranty

The warranty period for Miltex container gaskets is three years from the date of purchase.

Periodic Maintenance: Periodic Maintenance is suggested to ensure all parts are in good working order.

Repair or Replacement: Under this Warranty, and at its discretion, Miltex will repair or replace any parts of the product found to be defective in workmanship or material.

Miltex will determine, at its discretion, the cause and nature of any defect, the necessity and manner of repair or replacement, and all other matters pertaining to the conditions of the Miltex Container System.

Miltex, Inc. Returned Goods Policy applies to all returned items (refer to details in current price list).

Exclusions: This Warranty shall not apply to any condition(s) or damage resulting from negligence, misuse, improper handling, improper cleaning, improper opening techniques, or unauthorized repair work including but not limited to:

- Incidents of abuse such as denting of the Miltex Container due to dropping or other instances of mechanically applied pressure.
- Damage determined to be related to caustic or abrasive cleaning agents.
- Items modified by the customer.
- Items modified or customized by Miltex at the request of the customer.
- Damage from fire, flood and other occurrences not under the control of Miltex, Inc.



www.miltex.com
toll-free phone 866-854-8300
toll-free fax 866-854-8400
phone 717-840-9335
fax 717-840-9347

**Please call Customer
Service today for additional
information, an in-office
demonstration, to place an
order or to receive other
Miltex catalogs.**

**Illustrations and content
provide general description
only and may be subject
to change.**

Miltex, Miltex Redefining Excellence and EZ-Zyme are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation.
Silitone, Tape N'Tell, Tip-It, and Tuck-Away are trademarks of Miltex, Inc.
The Padgett logo is a registered trademark of Integra LifeSciences Corporation. Padgett is a trademark of Integra LifeSciences Corporation.
The Integra wave logo is a trademark of Integra LifeSciences Corporation.
Copyright © 2009 Integra LifeSciences Corporation

